Fall Facility Managers Meeting

November 1, 2016

Paul Umbeck, Director

Environment, Health & Safety Department Facilities Planning & Management University of Wisconsin-Madison

Agenda

- & Emergency Showers & Eyewash Stations
 - a Inspections & Maintenance
 - a Laboratories
 - A Hallway/Common Corridors
- & Laboratory's Role & Responsibility
- & EH&S' Role & Responsibility
- & Facility Manager's Role & Responsibility



October is Radiation Safety Refresher Training Month!

Click here for more information

October Radiation Refresher Training Month

External Links

HP Assist Login

Animal Risk Questionnaire

Pathogen Safety Data Sheets

About EH&S

Safety is a daily practice at UW-Madison. EH&S envisions a campus where safety awareness and compliance are integrated into everyday operations, making this University a model for safe research and education.

>> Learn More

POPULAR SEARCHES

Services

Abatement/Asbestos

Biosafety Cabinet Certification

Bio Protocol-ARROW

CORD

Dosimetry

Fire Extinguishers

Lab Remodeling Assessment

Lab Inspections/Outreach

Lake Conditions

Personal Protective Equip

Safety Data Sheets

Shower and Eyewash Stations

Guidance/Policy

Animal Contact

Research Bloodborne Pathogens

Chemical Hygiene Plan

Combustible Storage

DEA Controlled Substances

Facility Exit and Evacuation

Hazard Communication

Hazmat Shipping

Lab Safety Guide

Pesticides

Sharps Policy

Spill Information

QUICK LINKS



Emergencies

Protocol & contact info for emergency situations



Training

Descriptions & schedules for our training courses



Disposal Services

Disposal of chemicals, animal tissue & sharps



Forms/Manuals/Signs

Pull any form, manual or sign from our database



Contact List

Look here for staff member contact information

NEWS

Chemical Safety Training Course Now Offered at Multiple Locations

Our Chemical Safety Training course is now offered at more locations on campus, including Memorial Union and Union South. Click here to register!

Changes to BSC Billing Process

The BSC billing process has changed as of January 2015. Click here to learn more!

>> More News

What You Will Find on EH&S Website

- & Safety Overview
- & EH&S Annual Inspection
- k Keep the Area Clear
- & Self Inspection Guide
- & How to Report a Problem
- & Shower Code

Some Facts

- & UW-Madison has over 2,200 emergency showers/eyewash stations
- & EH&S/ Fire & Life Safety conducts an annual inspection and test of all showers and eyewash stations
- k The inspection findings are logged into our Fire & Life database
- Real There is no fee or chargeback to the department for this service
- & Plumbing related issues are submitted to the Plumbing Shop

Safety Overview

- & Laboratory eyewash stations should be checked weekly to ensure water flow and water quality are maintained
- & Areas around showers and eyewashes must be clear of clutter
- Eyewash nozzles should be protected by plastic caps (protects from dust and chemicals)
- NOT a substitute for PPE (safety glasses goggles, gloves, etc.)
- & Portable squeeze bottle are NOT acceptable as eyewash stations

Laboratory Role & Responsibility

- & Occupants of lab or adjacent room(s) should perform weekly testing of eyewashes in laboratories AND,
- & Should perform weekly testing of eyewash stations in hallways and corridors

Water Flow Appearance at an Eyewash Station



- when testing an eyewash station, water flow from each nozzle should be approximately six inches high and cross in the center
- Both streams of water should be equal

Which Eyewash has Good Flow?









The Winner is...



- The flows are equal and cross near the apex at the center of the eyewash
- The flows are tall enough that the user would not get their eyes too close to the nozzles

Facility Manager's Role & Responsibility

- k Track complaints or problems with eyewashes and showers
- & Submit work orders for any repairs that may be needed
- Educate, educate and educate staff about the importance of maintaining eyewash stations and emergency showers
- & Call EH&S Fire & Life Safety if you have questions